



# Product data sheet Granusan

V 02 / April 16, 2019

# Wetrok Granusan

Sanitary maintenance cleaner



#### Recommended dosage (in cold water):



**Maintenance cleaning**  
1 stick for 4 l of water



**Foam cleaning by hand**  
1 stick for 1 Wetrok foam bottle

#### Benefits:

Innovative acidic granulate cleaner for limescale removal and general cleaning. Suitable for medium hard water. More savings in storage and transport costs plus 100% dosing accuracy.

#### Application:

For maintenance cleaning of washable and acid-resistant surfaces (sanitary porcelain, fittings, ceramic tiles, bathrooms, showers, washrooms and WCs, etc.).


#### Please note:

Add the cleaning granulate to a filled foam bottle, close and shake the foam bottle vigorously for 30 seconds to completely dissolve the granulate. Wetrok Granusan is not suitable for acid-sensitive surfaces (limestone, marble, clinker, aluminium, etc.). The exterior of the sanitary fitting must not be warm when cleaning.

#### Properties:

Pre-portioned cleaning granulate  
Dissolves very quickly  
Highly acidic (in application solution)  
Scented  
High-lather  
Dyed red

#### pH value of the application solution:



#### Important ingredients:

Surfactants, acids, fragrance, colouring

#### Storage / shelf-life:

Store in original container in a cool, well-ventilated location. Shelf-life of the sticks (un-opened): 48 months. Use mixture within 7 days.

#### Recommended application methods:

##### Damp wiping (by hand)

Soak your cloth with cleaning solution and wring it out well. Wipe the entire surface with the damp cloth.

##### Damp wiping (with foam, by hand)

Apply the cleaning foam in an "X" form on a folded cloth, then damp wipe the surface (use the fold technique).