



# PRODUCT DATASHEET

Disinfection - Wipe-down surface disinfection

## Wetrok SintOx 15%

More eco-friendly disinfection application using hydrogen peroxide - effective against bacteria, yeast fungi, moulds, enveloped and non-enveloped viruses





# Wetrok SintOx 15%



**pH value (Concentrate):** 0.7–1.7

**Important ingredients:**

Hydrogen peroxide 7.5 g/100 g

## APPLICATION

- Surface disinfectant based on hydrogen peroxide/silver compounds
- Applications:  
Healthcare, veterinary medicine, food industry and feedstuffs sector
- Suitable for disinfection of drinking water, as protective agent in cooling and processing systems (CH only)

## PROPERTIES

- Bactericidal, yeasticidal
- Fungicidal (yeast fungi and moulds)
- Virucidal (polio, adeno, norovirus)
- Against enveloped and non-enveloped viruses
- Does not need to be rinsed off

## GOOD TO KNOW

- Fully effective against bacteria, viruses yeast fungi and moulds after 30 minutes

## PLEASE NOTE

- Brass may discolour
- Self-shining dispersions may be affected by regular application
- When using on paintwork or facing materials, it is advisable to test the product in a concealed spot before use

## RECOMMENDED APPLICATION METHOD

### Surface disinfection

Disinfect the cleaned, dry surface using a pre-wetted or lightly wrung-out mop, applying light pressure. Apply a sufficient amount of disinfectant. The entire surface must be wetted. Observe the contact time. Wear gloves!

## DOSAGE AND APPLICATION TIME

| Effectiveness             | Testing methods            | Dosage (%) | Contact time |
|---------------------------|----------------------------|------------|--------------|
| Bactericidal, yeasticidal | EN 13697, dirty conditions | 40         | 30 min       |
|                           |                            | 66         | 15 min       |
| Fungicidal                | EN 1650, clean conditions  | 50         | 15 min       |
| Virucidal                 | EN 14476                   | 66         | 30 min       |
| Norovirus                 | EN 14476                   | 66         | 15 min       |
|                           | EN 16777                   | 66         | 30 min       |
| Adenovirus                | EN 14476                   | 66         | 15 min       |
|                           | EN 16777                   | 66         | 30 min       |

## STORAGE

Store in original container at 15 to 25°C and protected from sunlight. Keep container closed when not in use. Protect against frost.

*Please refer to the product labels for information about the shelf life before opening.*

Wetrok AG shall have no liability whatsoever for damage caused by incorrect use. Observe special cleaning and care instructions from the floor covering manufacturer.

METHODS  
finder

